




2023-2024 BROCHURE


Blue Light Technology



SafeBlue System

SafeBlue Electrophoresis System / SafeBlue Imager System







MBE-G-Y1



MBE-150-3









ME15-UV15-TH01







ME15-UV15-TH02

Cat. No.	MBE-150	MBE-150-PLUS
Light control	On / Off switch	
Output voltage	10-150V / 1V	
Output current	10-300mA / 1mA	
Output power	30W	
Output type	Constant voltage or constant current	
Timer	999 (min) with alarm / Continuous	
Operating temperature	Ambient to 40°C	
Electrophoresis tank& amber filter lid	N/A	Yes
Output power Output type Timer Operating temperature	30W Constant voltage or constant curre 999 (min) with alarm / Continuous Ambient to 40°C	

Cat. No.	MBE-IMG-CA
Camera type	24.2MP digital camera, Wi-Fi function
Aperture	f/3.5-6.3(IS STM Lens)
Shutter speed	1/4000 - 30sec. (available range varies by shooting mode)
Storage media	SD memory card, Wi-Fi to PC, Smartphone or Tablet
Darkroom dimension	(WxLxH)340x195x400mm

SafeBlue System

SafeBlue Electrophoresis System

SafeBlue Imager System

for all

Gel and Sample Preparation .

Electrophoresis Running

Sample Observation Image Capture

for all

Traditional Method

Gel Electrophoresis System + Power Supply + Shaker + UV Transilluminator + Gel Documentation System

4 in 1 feature allows you to run as you see in realtime motion (run, view, capture and cut gel) with simple and easy set up on the image capture/ analysis software

MBE-150-PLUS Innovative Electrophoresis System

- 1. Early detection of running error (run as you see)
- 2. Prevention for breaking gel- no more gel transferring (no more gloves!!)
 - 3. Time saving gel staining & de-staining
 - 4. Compact in size (no large lab area required)

Life Sciences Research

ight Illuminator Blue

BluView Transilluminator, MBE-200A

- Thinner and more lightweight body
- Easy to carry
- Aluminum alloy casing design
- Low heat dispersion
- Energy saving product
- 470nm harmless blue light for direct human contact

Cat. No.	MBE-200A	
Dimension(mm)(WxLxH)	Approx. 200x200x13.9	
Viewing Area(mm)(WxL)	153x153	
Blue Light Source	15W	
Illuminator Base Design	Flat Bed	
Blue Light Wavelength	470nm	
Automatic Shutdown	Approx. 6 min	
Material	Aluminum Alloy	
Power	DC 12V, 2A	
Weight	Approx. 760g	











 \in

BluView Transilluminator, MBE- 00

- Portable 470nm blue light
- Slim design fits in most of electrophoresis tanks
- Real time observation
- Efficient and early stage mistake detection
- · Ultra high light uniformity
- Fit in a majority of different brands of mini and midi size...

Cat. No.	MBE-300	
Dimension(mm)(WxLxH)	Approx. 86x170x25	
Viewing Area(mm)(WxL)	Approx. 112x74.6	
Blue Light Source	20W	
Illuminator Base Design	Flat Bed	
Blue Light Wavelength	470nm	
Automatic Shutdown	Approx. 6 min	
Material	Aluminum Alloy	
Power	DC 12V, 2A	
Weight	Approx. 338g	



MBE-300 + ME 10 Tank horizontal electrophoresis system



MBE-300 + MJ-105A







* Technical specifications subject to change without notice.



www.majorsci.com info@majorsci.com



Taiwan Headquarters

Headquarters:No. 156, Sec. 1, Guoji Rd., Taoyuan Dist., Taoyuan City 33061, Taiwan T/+886-3-3762878 F/+886-3-3761310

Factory:

No.19, Ln. 207, Huakang St. Bade Dist., Taoyuan City 33464 Taiwan T/ +886-3-3623319 F/ +886-3-3623133 Email: info@majorsci.com

Shanghai Office

Room 612, International business exhibition center, 9300 Hunan Road, Pudong, Shanghai, China National toll-free No.:400-823-9177 T/ +86-21-50795277 F/ +86-21-50795277

US Office

19959 Sea Gull Way Saratoga, CA 95070 U.S.A. T/ +1-408-366-9866 F/ +1-408-446-1107

India Office

D. No. 12-13-99, Satguru Apts. Extn. Street. No. 3, Lane No.1, Tarnaka Secunderabad – 500 017. India T/ +91-40-27001515 T/ +91-40-27001586

