



NEW

Chemiluminescence Imager

MS-ChemiDoc



Features

- ➔ 10" Full-Color Capacitive Touchscreen: Intuitive and Smooth Operation
- ➔ Dual Mode (Physical Buttons & Touch-Free Gesture Control): Wave to operate the tray and prevent contamination
- ➔ Smart Status Indicators: Monitor imaging progress from a distance via ambient breathing lights
- ➔ Chemiluminescence Detection: Western Blot, Southern Blot, Northern Blot and other blotting experiments using substrates such as HRP, AP, etc.

Specification

Cat. No.	MS-ChemiDoc
Resolution	9 Megapixels
Quantum Efficiency (QE)	92% @ 540nm
Maximum Image Area	145 x 145mm
Pixel Size	3.76 x 3.76 μm
Readout Noise	0.9 e ⁻
Cooling Temperature	-25°C @ Ambient 22°C
Chemiluminescence Sample Tray	Single-layer
Image Storage	Internal Hard Drive, USB



www.majorsci.com
info@majorsci.com

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd.,
Taoyuan Dist., Taoyuan City
330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Building #C1, 880,
JiangYang Nan Road, Baoshan,
Shanghai, china 200439
T/ 13585666474

India Office

D.No.7-143, 2nd Floor,
St.No.2, Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586

* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice