



3D Rotating Mixer with Adjustable Speed

MS-3D-RM

Specification

Model	MS-3D-RM
Voltage	100~240V, 50/60Hz
Power	4W
Setting time	1-999min/0-999s
Speed	10-40rpm
Operation angle	20-99°
Working Mode	Continuous/Timed
Movement Mode	Rotation, Shaking, Angle
Application tube	50mL (single piece) 4 wells 15mL (single piece) 13 wells 1.5mL (single piece) 24 wells
Machine size	Approx. 240*140*210mm
Weight	Approx. 850g
Certification	CE



Ordering Info

Cat. No.	Product Description
MS-3D-RM-UK	3D Rotating Mixer with Adjustable Speed, UK
MS-3D-RM-EU	3D Rotating Mixer with Adjustable Speed, EU
MS-3D-RM-AU	3D Rotating Mixer with Adjustable Speed, AU
MS-3D-RM-US	3D Rotating Mixer with Adjustable Speed, US

* All images are for reference only, actual products might differ from the pictures above.
 * Technical specifications subject to change without notice.

Opional Accessory

Cat. No.	Product Description
MS-3D-RM-001	O-RING, OD14xID7.8mm, MS-3D-RM

Innovative Life Sciences Tools

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
 Taoyuan City 330041, Taiwan
 T/+886-3-3762878
 F/+886-3-3761310

Shanghai Office

Room 612, International business exhibition center,
 9300 Hunan Road, Pudong, Shanghai, China
 National toll-free No.:400-823-9177
 T/ +86-21-50795277
 F/ +86-21-50795277

US Office

19959 Sea Gull Way
 Saratoga, CA 95070
 U.S.A.
 T/ +1-408-366-9866
 F/ +1-408-446-1107

India Office

D.No.7-143, 2nd Floor,
 St.No.2,Nagendra Nagar,
 Habsiguda, Hyderabad-500007.
 India
 T/ +91-40-27001515
 T/ +91-40-27001586



*Could be compatibled with MO-MINI.